



# 650 nm Liquid Stop Solution for TMB Microwell

---

## **Product Description:**

650 nm Liquid Stop Solution for TMB Microwell is a propriety formulation used to stop BioFX One Component TMB Microwell Substrate reactions in the blue range of 620 nm to 650 nm.

## **Hazard Identification:**

650 nm Liquid Stop Solution for TMB Microwell is an acidic formulation should be treated with the same precautions as any other acid. Please refer to the MSDS for additional information.

## **Product Stability, Storage and Specifications:**

650 nm Liquid Stop Solution for TMB Microwell is stable for a minimum of 3 years when stored at room temperature (25°C). Store at room temperature (25°C). Refrigerated temperatures will not harm the reagent.

## **Product Use:**

650 nm Liquid Stop Solution for TMB Microwell is used to stop BioFX One Component TMB Microwell Substrates without producing a spectral change from blue to yellow. Equal volumes of substrate and stop solution should be used.

These products are for research and manufacturing use only and are not intended for use in humans, therapeutic or diagnostic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.